

Bien Air VS. NSK

Bien Air Optima INT BASIC



Optima INT **BASIC**

3 x CA1541

@\$1,770.00

@\$939.00ea

Total Package price = \$4596.00
3 Highspeeds + 3 Lowspeeds!!!!

NSK Ti-Max



NSK Ti-Max
Retail=
\$3,190.00
(XL 95L included)
Package \$ =
\$1,495.00

XL 95L X25L
Retail= Retail =
\$1,350.00 \$1,300.00
QTY x 2
Package\$ = Package\$ =
\$1,049.00 \$879.00

Total package price = \$4472.00
3 Highspeeds, 1 Lowspeed

| FEATURE | NSK TI-MAX | OPTIMA INT BASIC |
|----------------------------|--------------|------------------|
| Anti-Heat mechanism | N | Y |
| Endo using 1:1 | N | N/A |
| Implant capability | N | N/A |
| No Air in Endo Mode | N | N/A |
| Motor Torque (NCM) | 3 | 6 |
| Integration into ANY unit | N | Y |
| Motor RPM Range | 2,000-40,000 | 100-40,000 |
| RPM With 1:5 | 10K-200K | 500-200,000 |
| Swivel tubing | N | Y |
| Dual Glass-rod optics | N | Y |
| Dual Motor Capability | N | N/A |
| Noise level in DB | 62 | 50 |
| High & Low with 1:5 | N | Y |
| Display vs. actual bur spd | 37,000 | 40,000 |
| Display vs. actual bur spd | 191,000 | 200,000 |