

O-So Pure Dental 2020 Retail Pricing

Item #	Short Description	2020 Retail Price	Long Description
Disinfection Products <i>Disinfecting Water Lines</i>			
ADG-3600M	ADG-3600-Replacement Module ONLY	\$1,995	For use of replacing the module only
ADG-3600B	Disinfecting Dental Unit <u>BOTTLE</u> Line System	\$2,336	Bottle fed dental unit line shocking system with chemical free, unlimited volume and any size practice shocking (115v outlet required)
ADG-3600D	Disinfecting <u>DIRECT</u> Dental Unit Line System	\$2,742	Direct fed dental unit line shocking system with chemical free, unlimited volume and any size practice shocking (110/15A outlet required)
ADG-DF	ADG-3600 replacement desiccant cartridge	\$62	ADG-3600 Replacement Desiccant Cartridge for ADG-3600M / ADG-3600B / ADG-3600D
ATMOS-O3	Atmosphere-O3 Air sanitizer and purifier	\$656	This air purifier system purifies the atmosphere of any 2000sqft area utilizing ozone technology. Requires only a single 110/15A outlet, compact, lightweight and hidden from sight.
ATMOS-RC	Atmosphere-O3 replacemnt cell	\$143	Replacement cell for the Atmosphere air purifier (replace annually)
Steralization Center (SC) Products <i>Steralization Center Purification Products</i>			
SC-100	Steralization Center-100 is an all inclusive water treatment center for the 8 chair practice	\$2,354	SC-100 produces 4 GPH low-mineral (HF) for single instrument washer, dental unit bottles & mineral-free (DI) water for autoclaves / sterilizers. Includes: 3 gal tank, DI dispense gun, chrome faucet, TDS monitor. Requires one 110V/15A outlet. Installs in sterilization center cabinet.
SC-1PM	SC-100 Replacement Cartridge Kit (#1,#2,#3 & #4)	\$301	SC-100 replacement preventive maintenance kit. Contains 1 ea. Sediment (#1), carbon (#2), HF (#3) & DI (#4) cartridge filters.
SC-1DI	SC-100 Replacement DI (#4) cartridges, Set of two	\$158	Set of two DI (#4) filters for SC-100. . Replace when indicated by DI water quality monitor.
SC-1IWK	SC-100 Instrument washer installation kit	\$357	Instrument washer kit for SC-100. Includes 6 gallon tank, washer/disinfector manifold & install fittings
SC-400	Steralization Center-400 is an all inclusive water treatment center for the +8 chair capacity Practice/Clinic	\$3,785	SC-400 produces 16 GPH/.25 GPM low-mineral (HF) for multiple instrument washers, disinfected for dental unit bottles & mineral-free (DI) water for autoclaves / sterilizers. Includes: 10 gal tank, DI dispense gun, chrome faucet, TDS monitor. Requires one 110V/15A outlet. Installs in cabinet or countertop.
SC-4PM	SC-400 Replacement Filter Kit (#1,#2,#3 & #4)	\$465	Replacement Filter Cartridge kit for SC-400. Contains 5 cartridges.
SC4DI	Single SC-400 Replacement DI (#4) Cartridge.	\$175	Single DI filter Replacement Cartridge for SC-400. Replace when indicated by DI water quality monitor.
MECHANICAL ROOM (MR) Products <i>Mechanical Room Installation Filtration Products</i>			
MR-10	Mechanical Room Point-of-Entry Water Filtration & Manual ONLY Bypass	\$4,464	MR-10: Point-of-entry water filtration and UV disinfection system with auto filter flush with manual bypass. 10 GPM continuous flow, utilized with all plumbing. Requires: 110V/15A outlets. For dental unit waterline disinfection ADG-3600D required.
MR-10E	Mechanical Room Point-of-Entry Water Filtration & Switched Solenoid	\$4,872	MR-10E: Point-of-entry remote kill switch water filtration and UV disinfection system with auto filter flush, manual bypass, remote water control valve (24VAC or 110VAC). 10 GPM continuous flow, can be utilized with all plumbing. (requires 110V/15A outlets). For dental unit waterline disinfection ADG-3600D required

O-So Pure Dental 2020 Retail Pricing

<u>Item #</u>	<u>Short Description</u>	<u>2020 Retail Price</u>	<u>Long Description</u>
MR1PM1	Replacement Kit: MR-10 & MR-10E.	\$587	Complete Replacement Kit for MR-10 & MR-10E.
MR1PM2	High Output Replacement Kit: MR-10 & MR-10E	\$601	Optional High Output Replacement kit for MR-10 & MR-10E
MR1UVL	UV lamp for MR-10 & MR-10E; replace annually	\$382	UV Replacement Only for MR-10 & MR-10E; replace annually
MR4PM	Standard Replacement Filter kit for MR-400.	\$465	Standard Replacement kit for MR-400. Contains 4 cartridges,UV bulb.
MR4DI	Single MR-400 Replacement DI Cartridge	\$145	One DI filter for MR-400. Replace when indicated by monitor
Mineral Mizer			
DWS-RPKS	Annual service kit for DWS-SADI.	\$296	Annual service kit for DWS-SADI.
DWS-RPKSD	Annual service kit for DWS-SADI-C or DWS-SADI	\$325	Annual service kit for DWS-SADI-C or DWS-SADI
DWS-LRDC2	Set of two 10" DI filters for DWS-SADI, DWS-DDI, DWS-HERO-RODI, DWS-HERO-RODIUV.	\$249	Set of two 10" DI filters for DWS-SADI, DWS-DDI, DWS-HERO-RODI, DWS-HERO-RODIUV.
DWS-ROM-75HE	Membrane for DWS-SADI	\$156	RO Membrane100. For DWS-SADI, DWS-008.
DWS-ROM-100	Membrane for DWS-SADI, 008	\$199	RO Membrane100. For DWS-SADI, DWS-008.
DWS-UV9	UV lamp for DWS-001, DWS-004, DWS-005, DWS-008, DWS-HERO-RODIUV.	\$189	UV lamp for DWS-001, DWS-004, DWS-005, DWS-008, DWS-HERO-RODIUV.
DOMINATOR FAMILY			
DWS-RPKC	Dominator annual service kit	\$310	Dominator annual service kit
DWS-ROM-100	Dominator RO membrane (replace every 3 years)	\$199	Dominator RO membrane (replace every 3 years)
DWS-LRDC2	Mini-Mizer annual / interim DI filter service kit	\$249	Mini-Mizer annual / interim DI filter service kit
HERO FAMILY			
DWS-RPKH1	HERO Platinum/Complete annual service kit	\$505	HERO Platinum/Complete annual service kit
DWS-RPKH2	HERO Gold annual service kit	\$398	HERO Gold annual service kit
DWS-RPKH4	Annual service kit for DWS-HERO2-RODI-UV.	\$763	HERO Gold annual service kit
DWS-RPKH5	Annual kit for DWS-HERO2-RODI.	\$687	HERO Gold annual service kit
DWS-10DIC	Replacement DWS-HERO - DI	\$131	HERO Gold annual service kit
DWS-ROM-100H	RO Membrane, 100 GPD, DWS-HERO	\$270	HERO Gold annual service kit
PREDATOR FAMILY			
DWS-RPKD	Predator annual service kit	\$669	
DWS-2505	Standard 1st filter in MR10, DWS-POE, DWS-EF3.	\$153	Standard 1st filter in MR10, DWS-POE, DWS-EF3.
DWS-2055	Standard 2nd filter in MR10, DWS-POE, DWS-EF3.	\$280	Standard 2nd filter in MR10, DWS-POE, DWS-EF3.
DWS-2025	Hi-flow 2nd filter in MR10, DWS-POE, DWS-EF3. Use for high flow needs	\$234	Hi-flow 2nd filter in MR10, DWS-POE, DWS-EF3. Use for high flow needs.
BARRICADE FAMILY			
DWS-RPKL	Barricade annual service kit	\$305	Barricade annual service kit
TANK OPTIONS			
DWS-DIDG	DI water dispensing gun for sterilizer.	\$398	DI water dispensing gun & wall bracket for sterilizer.
DWS-ROSTB-14	10 gallon water storage tank	\$382	10 gallon water storage tank
DWS-ROSTB-14C	10 gallon high-flow water storage tank	\$473	10 gallon high-flow water storage tank
DWS-ROSTW-4	3.2 gallon water storage tank	\$164	3.2 gallon water storage tank
DWS-SFW	Chrome service faucet for dispensing self-contained bottle & drinking water.	\$91	Chrome service faucet (no air-gap) for dispensing self-contained bottle & drinking water.
DWS-ROSTADEC	3 gallon composite water storage tank for use in ADEC cabinets.	\$122	3 gallon pressurized rubber bladder water storage tank for use in ADEC cabinets.