

iOptima INT

6X WINNERS BEST OF CLASS



Whether you are planning to upgrade your dental practice with latest technology or you are just about to start your professional career: iOptima^{INT} is the All-in-One solution offering you the full range of treatment options seamlessly integrated into your dental chair. No need for additional equipment anymore: iOptima^{INT} lets you perform your restorative procedures, endo treatments (incl. reciprocating mode) and surgical procedures incl. implant placement. A Piezo scaling function is included as well.

- All-in-one device with special prep mode and endo mode with torque limit control
 - “No Air” endo Mode with Solenoid Valve
 - Complete library of endo files with presets
- Unlimited memory presets available in Prep and Endo Mode
- 3.5 Ncm MX2 Micromotor (the most powerful, stable motor ever measured) with ISO 3694 coupling
- Up to 4 motor option, with piezo, surgical motor and sterile irrigation pump
- Smartlogic sensorless, 3-phase voltage control allows for maximum torque at 100 RPM – eliminating the need for reduction attachments (like 7LP handle and 3:1 head) for restorative work
- Perfectly pairs with all attachments
- Apple display – the most intuitive display in the world
- Apple iPod or iPad mini – the largest display available in the industry
- Gorilla Glass® display
- System displays bur speed, live read out of real-time Torque, when selected gear ratio of handpiece is set
- Complete internal integration completely removes clutter and, “Breaking points”



ELECTROmatic™ Premium Electricis

Electricis Made Easy



Recipient of Dental Product Shopper Best Product Award.

Rated: 4.5/5

The ultimate all-in-one device for intuitive operation and integrated endo function. KaVo ELECTROmatic™ Premium, available with one or two motors for easy prep and integrated endo function. Available with one or two motors for easy prep and integrated endo function. The preprogrammed file systems make it easy to plan and customize your endodontic sequence.

- All-in-one device with special prep mode and endo mode with torque limit control
- Up to 5 speed memory presets available in prep mode
- Powerful KL703 LED Motor with ISO 3694 coupling
- Dual motor option for easy switch between high and low-speed functions
- SMARTdrive™ sensorless control allows a high torque at low speed ranges
- Pairs perfectly with KaVo MASTERmatic attachments
- Very easy to adjust speed with a slider on display
- Large 4.3" TFT color touch screen display
- Intuitive programming adjustment
- Scratch-resistant Dragontrail® glass display, which is six times harder than standard glass
- Smooth surface and easy to clean
- System displays bur speed when selected gear ratio of handpiece is set
- Remote mount display for best ergonomic viewing positions

Promo Comparison

Advanced Systems

Package W/ 2 high and Low



KAVO Part#	KAVO	KAVO \$	NSK Part#	NSK	NSK \$	MDW Part#	MDW	MDW \$	BA Part#	BA	BA\$
ElectroTorque Premium 1 motor		\$2,969.99	NLZ-E S120		\$3,379.00	E Motor		\$3,379.99	iOptima EXT		\$3,148.00
MasterMatic LUX M25L	 <small>MASTERmatic LUX M25 L High-speed Attachment™</small>	\$1,479.99	Ti-MAX Z95L		FREE	E Pro 1:5		FREE	EVO.15 1:5		\$1,714.00
MasterMatic LUX M25L	 <small>MASTERmatic LUX M25 L High-speed Attachment™</small>	FREE	Ti-MAX Z25L		FREE	E Pro 1:5		\$1,449.99	EVO.15 1:5 L		FREE
MASTERmatic LUX M07 L		\$1,069.00	NLZ Endo		FREE	E Pro 1:1		\$1,388.99	CA 1:1		\$1,429.00
INTRA LUX L80 1:1 Latch Low Sp		\$529.99	Ti-MAX Z95L		\$1,499.99	E Plus 6:1 Endo		\$1,159.99			
INTRA LUX L62 10:1 Endo		\$549.00									
Total		\$6,596.98			\$4,878.99			\$7,378.96			\$6,290.00

S W I S S  M A D E